



Detection of Bacterial Canker Pathogen in Tomato Seeds and Latent Infections of Seedlings

Dr. Hasan Bolkan
Former Director of Campbell
Vegetable R&D Davis California, USA

1



2




4



Agricultural Consulting



11



Seed Assay

12

Seed Assay

NO procedure is 100% correct

Use **Negative** and **Positive** samples for **most correct results**

Agricultural Consulting 13

Crop Losses due to Seed-borne Diseases

Immediate losses

- Reducing germination, poor seedling vigor,
- Abnormal seedlings, first crop losses from infected seedlings

Long term losses

- Pathogen surviving in soil
- Future inoculum source to new crops

Agricultural Consulting 14

Seed Health Testing: A Balance Act

– Management tool for the control of seed-borne/seed-transmitted diseases

- Increased product liability
- Competitive pressure within the industry
- Seed health becomes an important quality trait

Agricultural Consulting 15

Variation in Test Results

Sampling Error

Assay Sensitivity

Misdiagnosis & Saprophytes

Agricultural Consulting 16

Challenges in Detecting Seed-borne Pathogens



- Assay specificity/identification** • Need high level of confidence
- Contamination** • Target pathogen vs saprophytes
- Inoculum density** • Low level of the pathogen in seeds
- Sampling Error** • Too small of a sample/not representative sample

Agricultural Consulting 17

Challenges in detecting Seed-borne pathogens

1. A specific **method**/condition may be ideal for a given sample but not for another
2. Percent recovery is dependent on the **inoculum density**
3. Naturally infected seeds must be used for final **method evaluation**




Agricultural Consulting 18



Detection limit of the assay

- If detection limit is 1 contaminated seed in 10,000 healthy seeds: You need 10,000 seed replications (sample size is 30,000 seeds 3 X 10,000 seeds reps)

Agricultural Consulting 21



Conventional methods

Conventional

- Grow-out test (seedling test)
- Inoculating tomato plants with ground seed filtrates

Conventional

- Selective/differential media test
- General agar media test

Conventional

- Immunofluorescence staining (most commonly used)

Agricultural Consulting 25

Biotechnology based tests

Immunoassay based assays

- ELISA
- Poly/Monoclonal antibodies

Nucleic acid based assay

- DNA probes
- Classic PCR and Real-time PCR/Bio-PCR

Bacteriophage

- Phage Plaque Counting

Agricultural Consulting 26

PCR-Based Assay

Classic PCR: Culturing the pathogen from the infected seeds, extracting the DNA from the culture and subjecting it to PCR

Bio-PCR/Direct-PCR: Extraction of DNA from infected seeds and detection of the pathogen from the total population of bacteria on semi-selective media

Agricultural Consulting 27

Detection of Plant-Pathogenic Bacteria in Seed and Other Planting Material

Second Edition 2018

Edited by
M'Barek Fatmi,
Ron R. Walcott,
and Norman W. Schaad

APS PRESS
Agricultural Consulting 28

CHAPTER 17

Detection of *Clavibacter michiganensis* subsp. *michiganensis* in Tomato Seeds

M'Barek Fatmi, Hasan Bolkan, and Norman W. Schaad

Agricultural Consulting 29

U.S. National Seed Health System and ISHI Approved Standard Method

Incubation at 4C for 14 h

↓

Run PCR and/or pathogenicity test

Extraction using Stomacher

↙ ↘

Transfer suspected bacteria colonies on non-selective media


Dilution plating on semi-selective media

↙ ↘

Time needed 10-12 Days

Agricultural Consulting 30

The BIO/Real-Time PCR Protocol



- SCM**
 - Semi-selective Agar Medium
 - Ref: Fatmi & Schaad. 1988. *Phytopathology* 78: 121-126
- MSCM**
 - Modified Semi-selective Medium
 - Ref: Waters & Bolkan. 1992. *Phytopathology* 82: 1072
- SCM93**
 - Modified Semi-selective Medium 93
 - Ref: Fatmi & Atik. 2000. *Quatrieme congrès de l'AMPP, Rabat, Morocco*

Sterilized bag




The Stomacher



To Spike or not to Spike

Spiking of seed extracts with a known CMM culture to improve reliability of Cmm test



Courtesy of: Dr Harrie Koenraad

CMM colonies on different media




Colonial morphology of Cmm on different media. A: Typical colony on SCM after 9 days of incubation; B: Typical colony on mSCM after 5 days of incubation; C: Typical colony on NBY after 4 days of incubation.


Problems with PCR

- High sensitivity (contamination issues can give false positives)
- Lab must be very clean
- Operator very careful
- Chemicals of high purity (especially the tag polymerase)
- DNA, Primer & Mg concentrations
- Detection levels may not correlate to disease development

Pathogenicity Test



Conclusions



Detection of Latent Infections of Seedlings



LATEN INFECTION

- A stage in which a host is infected with a pathogen but does not show symptoms
- The interval from infection to display of symptoms (*incubation period*)

Detection Methods

Bioassay	<ul style="list-style-type: none">• Moist Blotters• Selective media
Histological	<ul style="list-style-type: none">• Microscopic examination
Serological & Molecular Techniques	<ul style="list-style-type: none">• ELISA• Real-time PCR/Classic PCR

Innovative Detection Methods

Spectroscopy-based methodology

- Analysis of volatile compounds as biomarkers

Remote sensing technologies

- Use of high-resolution portable spectral sensors to detect diseased leaves of tomato

Biosensors & Biophotonics

- devices that convert a biological response into an electrical signal
- detection of photons (quantum units of light)

Confirmation of CMM in the field by ImmunoStrips

Use of immunoStrips specific to CMM is a quick way to confirm the presence of bacterial canker disease in the field. However, be aware of false positives; confirmation of the disease by a laboratory may be needed. Sample preparation A, B, C



A - ImmunoStrips


B - Sample bag with buffer

C - Add suspected tissue sample and grind

(Photos copied from Agriplex brochure)

ImmunoStrip Assay

ImmunoStrip immersed in suspected ground tissue in buffer (A), positive results (B).



Latent infection



THANK YOU